

The orthopaedic knee walker Stable, adjust-NITRO GLIDE able handlebar Injuries, strain or fractures in the parking brak area of the lower thigh, ankle or metatarsal can significantly limit freedom of movement. In such situations our extremely agile **Nitro Glide** can offer an excellent solution. Designed as a knee walker, it enables you to move around quickly and safely, even when your legs are impaired. Corrosionresistant, robust frame and sporty

Smartsoft

tires for a

smooth ride

Removable

basket for

easy transport

of items

Holders for easy

transportation

of the crutches

After foot injuries, for example, with dislocations, sprains, fractures, or injuries in the leg, foot or ankle area.



INTENDED

Ergonomically-

shaped knee

As relief for temporary restrictions on the foot or lower leg, such as inflammation or swelling.



During the rehabilitation phase, after operations on the leg, foot or ankle.

## **TECHNICAL SPECIFICATIONS**

Item number	700500100
HMV no.	22.50.02.0XXX (in application)
rec. height	1,40 - 2,00 m
Total length	78 cm
Total width	40 cm   45 cm   50 cm
Handle height	78 - 99 cm
Knee pad height	42,5 - 57,5 cm
Knee pad height x width	34 x 17,5 cm
Folded dimensions	78 x 40 x 42,5 cm
Weight	10 kg
Max. Load	158 kg

